

## **Ti:Sapphire Laser Operating Procedure**

### **Start-up:**

- Step 1: Turn on water chillers using switch on back or side of each unit
- Step 2: Turn on pump laser using switch on back
- Step 3: Wait for the pump laser to warm up and temperature to reach equilibrium, at least 30 minutes
- Step 4: Put on laser goggles and turn on laser interlock
- Step 5: Turn on pump laser by turning key and opening shutter
- Step 6: Give Ti:Sapphire time to reach thermal equilibrium, at least 10 minutes

### **Optimization and Locking:**

- Step 7: Optimize light from pump laser into Ti:S
  - Measure the power after the Ti:S with a power meter at the beam splitter
  - Adjust the x and y alignment at the beam alignment unit between the pump laser and the Ti:S (inside the red box) to maximize power
- Step 8: Check etalon error on oscilloscope
  - Select "Etalon Error" from monitor A dropdown box on computer
  - Run a continuous scan on the resonator
  - The signal should cross zero, if it does not, repeat Step 7
- Step 9: Apply etalon lock on computer
- Step 10: Check resonator error on oscilloscope
  - Select "Resonator Error" from dropdown box
  - Run a continuous scan on the resonator
  - Signal should cross zero and peaks should be evenly spaced, if peaks are not evenly spaced laser is not operating on a signal spatial mode, unlock and re-lock etalon
- Step 11: Apply reference cavity lock on computer
- Step 12: Optimize light from Ti:S into the ECD-X
  - Look at "ECD PD 1" on oscilloscope
  - Adjust the x and y alignment between Ti:S and ECD-X to maximize signal on oscilloscope
  - Look at "ECD PD 2" on oscilloscope and optimize, it should optimize at the same x and y alignment as PD 1

- Alternately look at laser output from the ECD-X and adjust  $x$  and  $y$  alignment to narrow beam in each direction
- Step 13: Check ECD error on oscilloscope A: Error should cross zero B: If error does not cross zero, open “ECD” box on the computer and adjust “PD 1 gain” and “PD 2 gain” until signal cross zero
- Step 14: Apply ECD lock on computer A: If ECD does not lock, try increasing “Aux PD gain” in ECD box B: If ECD still does not lock, unlock, re-optimize and re-lock etalon and reference cavity
- Step 15: Optimize light from ECD-X into ECD-X-Q
  - Switch to other control console on computer
  - Look at “ECD PD 1” on oscilloscope
  - Adjust the  $x$  and  $y$  alignment between ECD- X and ECD-X-Q to maximize signal on oscilloscope
  - Look at “ECD PD 2” on oscilloscope and optimize, it should optimize at the same  $x$  and  $y$  alignment as PD 1
  - Alternately, look at laser output from the ECD-X-Q and adjust  $x$  and  $y$  alignment to narrow beam in each direction
- Step 16: Apply ECD-X-Q lock

**Shut-off:**

- Step 17: Remove all locks
- Step 18: Turn pump laser key to standby and close shutter
  - If laser will be used soon leave in standby
- Step 19: Cool the pump laser for a couple of hours
- Step 20: Turn off pump laser
- Step 21: Turn off water chillers